

[Home](#) > [Protocols](#) > Protocol for LunaScript® RT SuperMix Kit (E3010)

Protocol for LunaScript® RT SuperMix Kit (E3010)

1. Mix components briefly and spin down if necessary.
2. Prepare cDNA synthesis reaction as described below:

COMPONENTS	20 µl REACTION	FINAL CONCENTRATION
LunaScript RT SuperMix (5X)	4 µl	1X
RNA Sample	variable	(up to 1 µg)*
Nuclease-free Water	to 20 µl	

For no-RT control reaction, mix the following components:

COMPONENTS	20 µl REACTION	FINAL CONCENTRATION
No-RT Control Mix (5X)	4 µl	1X
RNA Sample	variable	(up to 1 µg)*
Nuclease-free Water	to 20 µl	

*Up to 1 µg total RNA, 1 µg mRNA or 100 ng specific RNA can be used in a 20 µl reaction. However, the cDNA input for downstream qPCR detection should typically contain < 10⁹ copies of the target to ensure that quantitation remains linear. To accommodate larger amounts of input RNA (> 1 µg), the reaction should be scaled up to ensure optimum cDNA synthesis.

For no template controls, mix the following components:

COMPONENTS	20 µl REACTION	FINAL CONCENTRATION
LunaScript RT SuperMix (5X)	4 µl	1X
Nuclease-free Water	to 20 µl	

[Cookies Settings](#)

[Accept cookies](#)

Incubate reactions in a thermocycler with the following steps:

CYCLE STEP	TEMPERATURE	TIME	CYCLES
Primer Annealing	25°C	2 minutes	1
cDNA Synthesis	55°C	10 minutes	
Heat Inactivation	95°C	1 minute	

Links to this resource

Related Products: [LunaScript™ RT SuperMix Kit](#)

Cookies Settings

Accept cookies

We use cookies to understand how you use our site and to improve the overall user experience. This includes personalizing content and advertising. To learn more and manage cookies, please refer to our [Cookie Statement](#).